

Procarta® Transcription Factor Plex Kit- By Request

Product Insert for TF Panel 2

Cat. Nos. PC5900–PC5939

About the Procarta Transcription Factor Plex Kit

Procarta Transcription Factor assays combine two technologies: Panomics' patent-pending detection probe separation method and xMAP® (multi-analyte profiling beads) from the Luminex® Corp. Together, these two technologies enable the quantitative measurement of DNA binding activity of up to 43 transcription factors from as little as 250 ng of protein from a nuclear extract or 500 ng of protein from a whole cell lysate.

Panomics' By Request Procarta Transcription Factor Plex Kits are available in 1-plate 96-well formats that can have a range from 2-42 plex. Each kit contains all of the reagents and assay controls required to prepare and detect DNA binding activity of customer specified transcription factors from human, mouse, or rat nuclear extracts. The Procarta Transcription Factor Whole Cell Lysis Kits is sold separately.

Safety Warnings and Precautions

CAUTION All chemicals should be considered potentially hazardous. We recommend that this product and its components be handled by those trained in laboratory techniques and be used according to the principles of good laboratory practice. This kit contains small quantities of sodium azide. Sodium azide reacts with lead and copper plumbing to form explosive metal azides. When disposing, flush drains with a large volume of water to prevent azide accumulation. Observe all state and local regulations for disposal. For research use only. Not for use in diagnosis of disease in human or animals.

Bead Analyte Association

Included (√)	Bead	Analyte	Included (√)	Bead	Analyte
	7	HAS+HBS		36	USF-1
	8	ALF-1/TAL-1		37	XRE
	9	TREF-1/2		38	GFI-1
	10	EGR		41	NPAS2
	11	NFκB		42	PIT-1
	12	c-MYC		43	SRY
	13	NeuroD1		44	E47
	14	PDX-1		45	CCAAT
	15	TR		46	AP-4
	17	MRE		47	CDP
	18	ELF-1		51	XBP-1
	19	SRE		52	CEF-1
	20	H4TF		53	ANTIOXIDENT RE
	21	HiNF		54	COUP-TF
	25	KPF-1		55	EVI-1
	26	LF-A1		61	GAG
	27	LvF		63	IKAROS
	28	MTF		65	HBS/XBP
	32	PUR-1		66	SIE
	33	RB		73	HSF
	34	TFE-3		76	TR (DR-4)
	35	AP-1			

Contents and Storage

Box 1 of 3 Shipped on Blue Ice Store entire box at 4 °C	Quantity 1-Plate Kit	Storage
Capture Beads, premixed	5.0 mL	2–8 °C
Streptavidin-PE	11 mL	2–8 °C
Elution Buffer	7.0 mL	2–8 °C
Binding Buffer	25 mL	2–8 °C
Separation Buffer	200 mL	2–8 °C
Assay Buffer	5.0 mL	2–8 °C
Wash Buffer	110 mL	2–8 °C
Reading Buffer	15 mL	2–8 °C
PCR Plate	1 each	2–8 °C
Sample Collection Plate	1 each	2–8 °C
Utility Plate	1 each	2–8 °C
Separation Plate	1 each	2–8 °C
Filter Plate	1 each	2–8 °C
PCR Plate Seal	1 x 5 each	2–8 °C
Plate Seal	1 x 5 each	2–8 °C

Box 2 of 3- Shipped on dry ice Store at -20 °C	Quantity 1-Plate Kit	Storage
Detection Probe, premixed	1.1 mL	-20 °C
Buffer A	20 mL	-20 °C
Buffer B	12 mL	-20 °C
DTT, 100 mM	350 µL	-20 °C
Protease Inhibitor	350 µL	-20 °C
Phosphatase Inhibitor I	350 µL	-20 °C
Phosphatase Inhibitor II	350 µL	-20 °C

Box 3 of 3- Shipped on dry ice Store at -80 °C	Quantity 1-Plate Kit	Storage
Positive Control Nuclear Extract	2 x 40 µL	-80 °C
Negative Control Nuclear Extract	2 x 40 µL	-80 °C
Positive Control Whole Cell Lysate	2 x 40 µL	-80 °C
Negative Control Whole Cell Lysate	2 x 40 µL	-80 °C

Contacting Panomics

For technical questions, contact our technical support group by telephone at 1-877-726-6642 option 3 or by email at techsupport@panomics.com (US and Canada) techsupport_europe@panomics.com (Europe), or visit our website www.panomics.com for an updated list of FAQs and product support literature.

Procarta is a registered trademark of Panomics, Inc. xMAP and Luminex are registered trademarks of Luminex Corp

© 2007 Panomics, Inc. All rights reserved.

Panomics Inc., 6519 Dumbarton Circle, Fremont, CA 94555, Tel. (877) 726-6642