

QuantiGene Technical Note

Recommendations for Defining a QuantiGene Plex Panel

About this Technical Note

This note provides information to help you define a QuantiGene Plex panel for your particular needs. It includes:

- ◆ Specifications for QuantiGene Plex 2.0 Magnetic Beads Assays
- ◆ Recommendations for selecting housekeeping genes
- ◆ Requirements for ordering a 2.0 Plex Set panel
- ◆ Panomics' available housekeeping genes for use in QuantiGene Plex 2.0 assays

QuantiGene Plex 2.0 Specifications

Parameter	Range
Limit of Detection	$\leq 1,000$ – $2,000$ transcripts/assay well
Limit of Quantitation	$\leq 2,000$ – $4,000$ transcripts/assay well
Linear Dynamic Range	≥ 3 logs
Assay Precision (CV)	$\leq 15\%$ intra-assay, $\leq 20\%$ inter-assay
Compatible Sample Types	Cultured cells, bacteria, whole blood, PAXgene blood or dried blood spots, fresh/frozen tissues (animal or plant), formalin-fixed paraffin embedded (FFPE) samples, purified RNA
Assay Format	96-well plate ^a
Targets/well	3–36

a. It is possible to run partial plates.

Selecting Housekeeping Genes

Every Plex Set should include 3–5 housekeeping genes. The purpose is to enable the following:

- ◆ Flexibility in data analysis, in case a selected housekeeper for a specific project does not behave as expected.
- ◆ Ability to normalize the data using a geometric mean of several housekeeping gene values. This technique “smooths” any bias from a single housekeeping gene.
- ◆ Provide alternatives, in case the signal from one housekeeping gene is too low or saturating.
- ◆ Flexibility for use for different sample types (cultured cells, blood, tissues, RNA).

Housekeeping Gene Characteristics

Try to select housekeeping genes that are expressed:

- ◆ Consistently, under all experimental conditions being evaluated.
- ◆ At a level similar to your target genes so that a single sample input results in a linear response for all targets.

Essential Information to Include in a Plex Set Order

Before you place an order for a Plex Set, make sure you have the following information:

- ◆ Unique identifier for each target gene, preferably an accession number AND version OR gi number. If neither of these is available, submit the exact nucleotide sequence of the transcripts you wish to detect.
- ◆ 3–5 housekeeping RNAs in each panel (**highly recommended**).
- ◆ Available information on relative level of expression of target and housekeeping RNAs in your experimental system, especially RNAs of low abundance.
- ◆ If data from another gene expression platform is to be compared to QuantiGene Plex 2.0 results, provide sequences of primers or probes used in that platform. We will design the Probe Sets to these regions whenever possible.
- ◆ If the panel is for validation of RNAi knockdown, you may provide the region(s) targeted by the siRNA. We will build the Probe Sets centered on these regions whenever possible.
- ◆ Additional specificity requirements. For example:
 - You wish to distinguish RNAs from human cells transplanted into a mouse model. In this case, we will design probe sets to be specific to the human RNA of interest but not the mouse homolog, and it will not cross react with any other human or mouse RNAs.
 - You wish to distinguish splice variants. In this case, provide both the exact sequence you wish to detect and the sequences of variants you do not wish to detect.

Note Unless otherwise specified, probe sets are specific to the sequence submitted or referred to by accession.version or gi number and will not cross react with any other RNAs (unless otherwise noted on the product insert) in the genome of the organism from which that RNA sequence is derived.

Relative Expression of Some Housekeeping Genes

The following table provides general guidelines for the linear assay working range based on the number of cells/assay well.

Housekeeping Genes and Expression Levels		Guidelines for Cells/Assay Well
Relative High Expression	Common Housekeeping Gene	2.0 Mag Capture Beads (Magnetic)
High	ACTB, B2M, GAPDH, PPIA, RPL19, RPS23, RPL32, RPLR0, RPS3, RPS18, RPS20	50–5,000
Medium High	LDHA, PGK1, PPIB, RPL13A, UBC	100–30,000
Medium	HPRT, POLR2A, TFRC	400–60,000
Low	ATP6V1A, GUSB, HMBS, TBP, TNX2	2,000–80,000

Housekeeping Gene Recommendations Based on Sample Type The following table provides recommendations for housekeeping genes based on the type of sample you are examining.

Housekeeping Genes and Expression Levels		Sample Types ^a				
Relative Expression Level	Housekeeping Gene	Whole Blood or PAXGene Blood	Cultured Cells	Fresh/Frozen Tissues	Purified RNA	FFPE Tissues
		Maximal Sample Input	20,000 Cells/Well	Maximal Sample Input	250 ng Total RNA	Maximal Sample Input
High	ACTB, B2M, GAPDH, PPIA, RPL19, RPS23, RPL32, RPLR0, RPS3, RPS18, RPS20	+++	+	++	+	++
Medium High	LDHA, PGK1, PPIB, RPL13A, UBC	++++	+++	+++	+++	++++
Medium	HPRT, POLR2A, TFRC	++	++++	++++	++++	++++
Low	ATP6V1A, GUSB, HMBS, TBP, TNX2	+	+++	++	+++	+

a. + = not recommended, ++ = fair, +++ = good, ++++ = best

Common Housekeeping Genes The following table provides the full name and accession number of the housekeeping genes listed in the previous tables.

Symbol	Name	Human Accession Number	Mouse Accession Number	Rat Accession Number
ACTB	Actin, beta	NM_001101	NM_007393	NM_031144
ATP6V1A	ATPase, H+ transporting, lysosomal 70kDa, V1 subunit A	NM_001690	NM_007508	XM_001058034
B2M	Beta-2-microglobulin	NM_004048	NM_009735	NM_012512
GAPDH	Glyceraldehyde-3-phosphate dehydrogenase	NM_002046	NM_008084	NM_017008
GUSB	Glucuronidase, beta	NM_000181	NM_010368	NM_017015
HMBS	Hydroxymethylbilane synthase	NM_000190	NM_013551	NM_013168
HPRT1	Hypoxanthine phosphoribosyltransferase 1 (Lesch-Nyhan syndrome)	NM_000194	NM_013556	NM_012583
LDHA	Lactate dehydrogenase A	NM_005566	NM_010699	NM_017025

Symbol	Name	Human Accession Number	Mouse Accession Number	Rat Accession Number
PGK1	Phosphoglycerate kinase 1	NM_000291	NM_008828	NM_053291
POLR2A	Polymerase (RNA) II (DNA directed) polypeptide A 220kDa	NM_000937	NM_009089	XM_343922
PPIA	Peptidylprolyl isomerase A (cyclophilin A)	NM_021130	NM_008907	NM_017101
PPIB	Peptidylprolyl isomerase B (cyclophilin B)	NM_000942	NM_011149	NM_022536
RPL13A	Ribosomal protein L13A	NM_012423	NM_009438	NM_173340
RPL19	Ribosomal protein L19	NM_000981	XM_001476576	NM_031103
RPL32	Ribosomal protein L32	NM_000994	NM_171086	NM_013226
RPLP0//Arbp	Ribosomal protein, large, P0 (human) // acidic ribosomal protein P0 (mouse and rat)	NM_001002	NM_007475	NM_0022402
RPS3	Ribosomal protein S3	NM_001005	NM_012052	NM_001009239
RPS18	Ribosomal protein S18	NM_022551	NM_011296	NM_213557
RPS20	Ribosomal protein S20	NM_001023	NM_026147	NM_001007603
RPS23	Ribosomal protein S23	NM_001025	NM_024175	NM_078617
RPS29	Ribosomal protein S29	NM_001032	NM_009093	NM_012876
TBP	TATA box binding protein	NM_003194	NM_013684	NM_001004198
TFRC	Transferrin receptor (p90, CD71)	NM_003234	NM_011638	XM_340999
TXN2	Thioredoxin 2	NM_012473	NM_019913	NM_053331

Contacting Panomics For technical support, contact the appropriate resource provided below based on your geographical location. For an updated list of FAQs and product support literature, visit our website at www.panomics.com.

Location	Contact Information
North America	1.877.726.6642 option 3 techsupport@panomics.com
Europe	+39.02.95.360.250 techsupport_europe@panomics.com
Asia	+86.10.59208157 techsupport_asia@panomics.com

QuantiGene® is a registered trademark exclusively licensed to Panomics, Inc.

© 2008 Panomics, Inc. All rights reserved.

Panomics Inc., 6519 Dumbarton Circle, Fremont, CA 94555, Tel. (877) 726-6642

Rev. A 081202