

**FAST-1 EMAS Kit****Product Insert****Cat. No. AY1015P**

---

**Amount** 25  $\mu$ L (10 ng/ $\mu$ L)

---

---

**Lot No.** Specified on product label

---

---

**Storage** -20°C

---

---

**Componets** A Biotin labeled probe intended for use with Panomics EMSA Kits. Each Labeled TF probe is accompanied with 25  $\mu$ L of Cold Probe (330 ng/ $\mu$ L), its unlabeled form.

---

---

**Description** FOXH1: forkhead box H1

---

---

**Sequence (5'->3')** CGGATTGTGTATTGGCTGTAC

---

---

**Reference** Mol. Cell 2:121-127 (1998)

---

---

**Contacting  
Panomics**

For ordering information or technical support, contact the appropriate resource provided below according to your geographical location.

<b>Location</b>	<b>U.S. Corporate Headquarters</b>
Address	Panomics, Inc. 6519 Dumbarton Circle Fremont CA 94555 USA
Telephone	1.510.818.2600
FAX	1.510.818.2610
Email	info@panomics.com
Technical Support	<a href="mailto:techsupport@panomics.com">techsupport@panomics.com</a> or <a href="tel:1.877.726.6642">1.877.726.6642 option 3</a>
Ordering Information	orders@panomics.com

<b>Location</b>	<b>European Headquarters</b>
Address	Panomics Srl Via Sardegna 1 20060 Vignate Milano Italy
Telephone	+39.02.95.360.250
FAX	+39.02.95.360.992
Email	info_europe@panomics.com
Technical Support	<a href="mailto:techsupport_europe@panomics.com">techsupport_europe@panomics.com</a>
Ordering Information	order_europe@panomics.com

<b>Location</b>	<b>Asia Pacific Headquarters</b>
Address	Panomics, Inc. 16F Gemdale Plaza Tower A, No. 91 Jianguo Road Beijing 100022 P.R. China
Telephone	+86.10.59208157
FAX	+86.10.59208111
Email	info_asia@panomics.com
Technical Support	<a href="mailto:techsupport_asia@panomics.com">techsupport_asia@panomics.com</a>
Ordering Information	order_asia@panomics.com

For an updated list of FAQs and product support literature, visit our website at [www.panomics.com](http://www.panomics.com).

---