

DeliverX™ siRNA Transfection Kit 1.0 mL

Product Insert

Cat. No. DX0003

About DeliverX siRNA Transfection Kit

Panomics' DeliverX siRNA Transfection Kit is suitable for the transfection of most cell types. The kit contains the reagents required to efficiently transfect siRNAs into most cell types with minimal damage and good reproducibility when optimization guidelines provided in the User Manual are followed.

Safety Warnings and Precautions

All chemicals should be considered potentially hazardous. We recommend that this product and its components be handled by those trained in laboratory techniques and be used according to the principles of good laboratory practice.

For research use only. Not for use in diagnosis of disease in human or animals.

Contents and Storage

Component	Volume	Storage
DeliverX siRNA Transfection Reagent	1.0 mL	-20 °C
DeliverX siRNA Buffer-1	34 mL	15-30 °C
DeliverX siRNA Buffer-2	34 mL	15-30 °C

- ◆ We recommend dividing DeliverX siRNA Transfection Reagent into 50–100 µL aliquots to minimize freeze/thaw cycles.
- ◆ If stored properly, product shelf life is 6 months from the date received.

Contacting Panomics

For technical questions, contact our technical support group by telephone at 1-877-726-6642 option 3 or by email at techsupport@panomics.com (US and Canada) techsupport_europe@panomics.com (Europe), or visit our website www.panomics.com for an updated list of FAQs and product support literature.

DeliverX is a trademark of Panomics, Inc.

DeliverX and DeliverX Plus transfection reagents are manufactured and distributed by Panomics, under license from CNRS (France).

© 2006 Panomics, Inc. All rights reserved.

Panomics Inc., 6519 Dumbarton Circle, Fremont, CA 94555, Tel. (877) 726-6642