

## TIMP3 Gene Promoter Reporter Vector

### Product Insert

**Cat. No. LR1034**


**Shelf Life** One year from date of receipt under proper storage conditions.

**Storage** Store at -20°C

### Contents

Gene Promoter Reporter Vector	20ul of 500ng/ul; 10ug
Empty Control Reporter Vector	20ul of 500ng/ul; 10ug
Sequence information is available on our website	

Gene Promoter	TIMP3
Description	Designed to study regulation of promoter for Homo sapiens TIMP metalloproteinase inhibitor 3 (Sorsby fundus dystrophy, pseudoinflammatory) (TIMP3)
RefSeq ID (for downstream gene)	NM_000362



### Sequence (5'-3')

```

AAAAAGATGTGGAAACTGGTTCCCTTTCCTATGGCTAAAGTGCTCAAAGGGGAAAAAGGAGTTTCAAAAATGTTCTTGG
AAATACCATCTCTCACGAATTCTTCGGCCTCTGCTGTCCCAATGTCACTTGTCTGAGATGTAAACAGAGGAGTTCTGAGA
AAGAAGCTGAACCTGCATTCTCCCTGTTTCTATTTGTTCCAAACTTGTGGCATTCTAACAGGATGAAGCGGAAGAGAA
AGGAAAAGAGACAAAAGTGTAAGAAAGATGGAAGATCCCAGCTGCAAATGGCCATTTCAGTTAGATGGAACAGCTGCTGA
CGTTCAGGGAAATGCATGTCTCTCTTCAGATGGGAAGGAGCAGTGGAAAAGGGGTGACGAGTTCCTGGCTGGCCACCAATC
ATCCCATCTTTCTGTGCCGTTTCTCATCTGGAAAGTGGGAGTGATACTTGTGCTTGCTTTTCTACCCACAAAGATTAT
TGTGAGAGCTATAATACGGTGAGATACAGAATCCTGCTTTTAAAAATACAAAGCAGAATCAAGATGTCAATAATAAGGAT
AGTAATTGTGTAGTTATCTGCAATCATCTATTATAGCTAGTCGTAGGATCCTGGATCGTTCTCCTGGTTTTACTACA
GTTTTGGATCAGCTCACCCCAAATCCCTTGTGAAGGGTGGAGCCCTGTCAGCCATGGGCAGGGAACCACTTCTCCTTGG
CCTTTCTACTTTCTGTCCCTTCAAACATGCCAGGGTCTTTGCACTTGTGTTCCCCCTGCCTGGTACCTCTCCTCTGTGG
CTTGCCCCAGAGCTGATCCTTGTCTTTGCCACTTCTCAGCGAGGATGGCACTTCAGGGAGCCCTTCCCTTACTATCGCA
GAGAGAGCAGGCCCTCCCAAGTCATGTCCAACCCAGAACTCTGTTTTGTTTTCTCATAGCCCTAGCATCACAGAAAATC
ACCCTGTGCATTCATGGATGTCCACGGGGCAAGGGCTTT

```

---

**Contacting  
Panomics**

For ordering information or technical support, contact the appropriate resource provided below according to your geographical location.

<b>Location</b>	<b>U.S. Corporate Headquarters</b>
Address	Panomics, Inc. 6519 Dumbarton Circle Fremont CA 94555 USA
Telephone	1.510.818.2600
FAX	1.510.818.2610
Email	info@panomics.com
Technical Support	<a href="mailto:techsupport@panomics.com">techsupport@panomics.com</a> or <a href="tel:1.877.726.6642">1.877.726.6642 option 3</a>
Ordering Information	orders@panomics.com

<b>Location</b>	<b>European Headquarters</b>
Address	Panomics Srl Via Sardegna 1 20060 Vignate Milano Italy
Telephone	+39.02.95.360.250
FAX	+39.02.95.360.992
Email	info_europe@panomics.com
Technical Support	<a href="mailto:techsupport_europe@panomics.com">techsupport_europe@panomics.com</a>
Ordering Information	order_europe@panomics.com

<b>Location</b>	<b>Asia Pacific Headquarters</b>
Address	Panomics, Inc. 16F Gemdale Plaza Tower A, No. 91 Jianguo Road Beijing 100022 P.R. China
Telephone	+86.10.59208157
FAX	+86.10.59208111
Email	info_asia@panomics.com
Technical Support	<a href="mailto:techsupport_asia@panomics.com">techsupport_asia@panomics.com</a>
Ordering Information	order_asia@panomics.com

For an updated list of FAQs and product support literature, visit our website at [www.panomics.com](http://www.panomics.com).

---